

Candidate's Name:.....

Signature:.....

P525/3

CHEMISTRY

(Practical)

Paper 3

Jul/Aug.2023

3 ¼ hours

WWEK PRACTICALS

Uganda Advanced Certificate of Education

CHEMISTRY-PRE MOCK SET7

(PRACTICAL)

Paper 3

3 hours 15 minutes

INSTRUCTIONS TO CANDIDATES:

*Answer all questions. Use **blue** or **black** ball point pens. Any work done in pencil will **not** be marked except drawings.*

Record your answers on this question paper in the spaces provided.

Mathematical tables and silent non-programmable scientific electronic calculator may be used.

*Reference books (i.e. text books, booklets on qualitative analysis etc.) should **not** be used.*

*Candidates are not allowed to start working with the apparatus for the first **15 minutes**. This is to enable candidates read the question paper and make sure they have all the apparatus and chemicals that they may need.*

For Examiners' Use only

Q.1	Q.2	Q.3	TOTAL

1. You are provided with the following:

FA1, which is a solution containing 39.2 g per litre of ammonium ferrous sulphate $(\text{NH}_4)_2\text{SO}_4 \cdot \text{FeSO}_4 \cdot 6\text{H}_2\text{O}$.

FA2, which is a solution of potassium manganate (VII) of unknown concentration.

2M sulphuric acid

Solid **W**

You are required to determine the;

- (i) Concentration of **FA2** in moles per litre
- (ii) Stoichiometric ratio of reaction between Fe^{2+} ions in **FA1** and **W**.

Manganate (VII) ions reacts with iron (II) ions in acidic medium according to the following equation.



PROCEDURE:

PART I:

Pipette 25.0 cm^3 (or 20.0 cm^3) of **FA1** into a conical flask, then add 20 cm^3 sulphuric acid. Titrate the mixture with **FA2** from the burette until the end point. Repeat the titration until you obtain consistent results. Record your results in **Table 1**.

Results:

Volume of pipette used..... cm^3 (1/2 mark)

Table 1

Final burette reading (cm^3)			
Initial burette reading (cm^3)			
Volume of FA2 used (cm^3)			

(4 1/2 marks)

Titre values used to calculate average volume of **FA2** used

(1/2 mark)

.....cm³

Average volume of **FA2** used

(2 1/2 marks)

.....
.....
.....

Questions:

Calculate the concentration of manganate (VII) ions in **FA2** per litre of solution. (H= 1, N= 14, O= 16, S=32, Fe= 56)

(3 marks)

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....
.....

PART II:

Weigh accurately about 0.5g of **W** and place it in a beaker. Add to it about 50cm³ of distilled water and stir to dissolve. Transfer the contents of the beaker into a 250cm³ volumetric flask. Add exactly 150cm³ of **FA1** and top up with distilled water to the mark, heat the mixture for about 10 minutes. Shake and allow to stand for about 5 minutes. Label the resultant solution **FA3**.

Pipette 25.0cm³ (20.0cm³) of **FA3** into a conical flask followed by 10cm³ of 2M sulphuric acid and then titrate with **FA2** until the end point. Repeat the titration until you obtain consistent results. Record your results in **Table 2**.

Results:

Mass of **W** and the weighing bottle.....g (1/2 mark)

Mass of empty weighing bottle.....g (1/2 mark)

Mass of **W** used.....g (1/2 mark)

Volume of pipette used.....cm³ (1/2 mark)

Table 2

Final burette reading (cm ³)			
Initial burette reading (cm ³)			
Volume of FA2 used (cm ³)			

(4 1/2 marks)

Titre values used to calculate average volume of **FA2** used (1/2 mark)

.....cm³

Average volume of **FA2** used (2 1/2 marks)

.....

Questions:

(a) Calculate the number of moles of;

- (i) Excess iron (II) ions that reacted with manganate (VII) ions in **FA2**. (2 ½ marks)

.....

.....

.....

.....

.....

.....

.....

(ii) Excess iron (II) ions contained in 250cm³ in **FA3**. (2 marks)

.....

.....

.....

.....

.....

.....

(iii) Iron (II) ions that reacted with **W**. (2 ½ marks)

.....

.....

.....

.....

.....

.....

(b) Determine the reaction ratio between iron (II) ions and **W**. (Relative formula mass of **W** = 123) (2 ½ marks)

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

2. You are provided with substance **X** which contains two cations and two anions. You are required to identify the ions in **X**. Carry out the following tests and identify any gas(es) evolved. Record your observations and deductions in **Table 3**. (30 marks)

TESTS	OBSERVATIONS	DEDUCTIONS
(a) Heat X in a dry test tube gently then strongly and allow to cool.		
(b) To one spatula endful of X add concentrated sulphuric acid and warm.		

(c) Shake two spatula endfuls of X with 5cm ³ of water. Divide the resultant solution into three parts.		
(i) To the first part, add 3-4 drops of ethanol followed by 2-3 drops of concentrated sulphuric acid and heat.		
(ii) To the second part, add Iron (III) chloride solution and heat.		
(iii) To the third part add ammonia solution dropwise until in excess. Filter, keep both the filtrate and residue.		
(d) To the filtrate from (c iii), add dilute hydrochloric drop wise until the solution is just acidic. Divide the resultant solution into five parts.		
(i) To the first part of the acidic solution, add sodium hydroxide solution dropwise until in excess.		
(ii) To the second part of the acidic solution, add ammonia solution dropwise until in excess.		

<p>(iii) Use the third part of the acidic solution to carryout a test of your choice to confirm the first cation in X.</p> <p><u>Test:</u></p>		
<p>(iv) To the fourth part of the acidic solution, add 2-3 drops of lead (II) nitrate solution.</p>		
<p>(v) To the fifth part of the acidic solution, add 2-3 drops of barium chloride solution.</p>		
<p>(e) Wash the residue from (c iii) with ammonia solution once then dissolve it in dilute sulphuric acid. Divide the resultant solution into three parts.</p>		
<p>(i) To the first part add sodium hydroxide solution dropwise until in excess.</p>		
<p>(ii) To the second part add ammonia solution dropwise until in excess.</p>		

(iii) To the third part add litmus solution followed by ammonia solution.		
--	--	--

Cations in **X** are.....and.....

Anions in **X** are.....and.....

2. You are provided with substance **M** which is an organic compound. You are required to determine the nature of **M**. Carry out the following tests and identify any gas(es) evolved.

Record your observations and deductions in **Table 4**.

(20 marks)

Table 4.

TESTS	OBSERVATIONS	DEDUCTIONS
(a) Burn a small amount of M on a spatula end		
(b) To 2cm ³ of M , add 3cm ³ of water and shake. Divide the mixture into three parts.		
(i) To the first part of the solution add 2-3 drops of Iron (III) chloride solution.		

(ii) To the second part of the solution add 2-3 drops of sodium carbonate solution.		
(iii) To the third part of the solution add 2-3 drops of acidified potassium manganate (VII) solution and heat.		
(c) To 0.5cm ³ of M in a test tube, add 2-3 drops of 2,4- dinitrophenyl hydrazine solution (Brady's reagent)		
(d) To 1cm ³ of M , add 2cm ³ of ethanol followed 2-3 drops of concentrated sulphuric acid and heat. Pour the products in a beaker of cold water and stir.		
(e) To 1cm ³ of M , add 2-3 drops of Tollen's reagent and heat.		

Describe the nature of **M**.

.....

.....

.....

.....

END

CONFIDENTIAL

Each student is to be provided with;

120 cm³ of **FA1**

60cm³ of **FA2**

1.0 g of solid **W**

60cm³ 2M Sulphuric acid

Weighing scale

100 ml measuring cylinder

3 beakers (250ml each)

Volumetric flask (250ml)

Distilled water Pipettes (25.0cm³/20.0cm³)

2 conical flasks

Burette

Droppers

Source of heat

A spatula

3 g of substance **X**

6 test tubes

2 filter papers

Litmus papers

FA1 is made by dissolving 39.2g of Ammonium ferrous sulphate in 1 litre

FA2 is 0.02M Potassium manganate (VII) solution

Solid **W** is Potassium chlorate.

X is a mixture of Zinc acetate and Aluminium sulphate in a ratio 3:2

3g of substance **W**.

M is Methanoic acid/Formic acid.