

DSF UACE BIOLOGY PRACTICAL CONFIDENTIAL JUNE 2023

Provide each student with the following:

Specimen **Z** (Freshly killed, mature toad/ frog).

Solution **X** (Hydrogen peroxide 3%) 100 ml.

Solution **Y** (fresh Irish potato extract prepared as provided hereunder) - 20 ml.

Distilled water - 250 ml.

5 Test tubes of equal size.

2 Beakers 100 ml.

10 cm ruler.

Sticky label paper.

Measuring cylinder - 50-100 ml.

Sieve with fine pores/ cloth.

20 small discs from a filter paper, made using a clean hole punch.

Stopwatch or digital timing device.

Specimen **R₁** (Mature tuber of Irish potato with shoots).

Specimen **R₂** (Mature tuber of Irish potato without shoots).

Specimen **S** (A few mature leaves of Bryophyllum with young adventitious buds).

Specimen **T** (Aquatic offset e.g., *Eichhornia*/ *Pistia* OR terrestrial stolon e.g., *paspalum*/ spear grass) having at least two shoots: one new, another older).

Hand lens (magnifying glass).

Dissection tray or dissecting pan.

Dissection tools (scalpel, scissors, forceps, and probe)

Paper towels or cotton wool.

Pins or T-pins.

How to prepare solution Y for 5 students:

Peel fresh Irish potatoes and wash them thoroughly.

Use a blender to crush 50g of Irish potato mixed with 50 ml distilled water.

Filter the potato extract through a sieve with fine pores, into a measuring cylinder.

Add distilled water into the filtrate to make a total volume of 100 ml.