

Name.....

Signature .....

**P525/3**

**CHEMISTRY**

**(Practical)**

**Paper 3**

**Sep/Oct. 2023**

**3 ¼ HOURS**

**WWEK PRACTICALS**

**Uganda Advanced Certificate of Education**

**S6 CHEMISTRY**

**(PRACTICAL)**

**Paper 3**

**3 hours 15 minutes**

**INSTRUCTIONS TO CANDIDATES:**

Answer **all** questions. Use **blue** or **black** ink. Any work done in pencil will **not** be marked **except** drawings.

*Record your answers on this question paper in the spaces provided.*

*Mathematical tables and silent non-programmable scientific calculators may be used.*

*Reference books (i.e. text books, booklets on qualitative analysis, etc.) should **not** be used.*

*Candidates are **not** allowed to start working with the apparatus for the first **15 minutes**. This time is to enable candidates to read the question paper and make sure they have all the apparatus and chemicals that they may need.*

FOR EXAMINER'S USE ONLY			
Q.1	Q.2	Q.3	TOTAL

1. You are provided with the following:

**FA1**, which is a 0.1M hydrochloric acid

**FA2**, which is approximately a 0.1 M sodium hydroxide solution

**FA3**, which is 2M Hydrochloric acid

Solid **W**, which is impure metal carbonate.

You are required to standardize **FA2** and use it to determine the percentage of the impurity in **W**.

**PROCEDURE I:**

(a) Pipette 20.0cm<sup>3</sup> or 25.0cm<sup>3</sup> of **FA2** into a conical flask. Titrate with **FA1** using phenolphthalein as the indicator. Repeat the titration until you obtain consistent results. Record your results in table below.

**Table:**

Volume of pipette used .....cm<sup>3</sup>

Final burette reading (cm <sup>3</sup> )			
Initial burette reading (cm <sup>3</sup> )			
Volume of <b>FA2</b> used (cm <sup>3</sup> )			

Values used for calculating average of **FA1**

.....and .....cm<sup>3</sup>

Average volume of **FA1**

.....  
.....  
.....  
.....  
.....  
.....

**Question:**

Calculate the concentration in moles per litre of solution **FA2**.

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

(b) Weigh 1.5g of **W** and dissolve it in 20.0cm<sup>3</sup> of **FA3** in a beaker. Add 80cm<sup>3</sup> of distilled water and label this solution **FA4**.

Pipette 20.0cm<sup>3</sup> or 25.0cm<sup>3</sup> of **FA2** into a conical flask and titrate it with **FA4** using phenolphthalein as the indicator. Record your results in the table below.

**Results:**

Mass of container + **W** = .....g

Mass of empty container = .....g

Mass of **W** = .....g

**Table:**

Volume of pipette used .....cm<sup>3</sup>

Final burette reading (cm <sup>3</sup> )			
Initial burette reading (cm <sup>3</sup> )			
Volume of <b>FA2</b> used (cm <sup>3</sup> )			

Values used to calculate the average of **FA4**.

..... and .....cm<sup>3</sup>

Average volume of **FA4** used

.....

.....

.....

.....

**Questions:**

(a) Calculate the:

(i) Number of moles of hydrochloric acid that was in 100cm<sup>3</sup> of **FA4**.

.....

.....

.....

.....

.....

.....

(ii) Number of moles of the metal carbonate that reacted. (the metal carbonate reacts with hydrochloric acid in the ratio 1:2)

.....

.....

.....

.....

(b) Determine the percentage of the impurity in **W**. (Relative formula mass of the metal carbonate is 84)

.....

.....

.....

.....

.....

.....

2. You are provided with a substance **P** containing two cations and two anions. You are required to carry out tests on substance **P** to identify the cations and anions in it. Identify any gas(es) evolved. Record your observations and deductions in Table 3 below. (30 marks)

**Table 3**

TESTS	OBSERVATIONS	DEDUCTIONS
(a) Shake <b>two</b> spatula endfuls of <b>P</b> with 5cm <sup>3</sup> of water. Filter, keep both the <b>filtrate</b> and <b>residue</b> . Divide the filtrate into <b>four</b> parts.		
(i) To the <b>first</b> part of the filtrate, add 2-3 drops of dilute sulphuric acid.		
(ii) To the <b>second</b> part of the filtrate, add 2-3 drops of lead (II) nitrate solution followed by sodium hydroxide solution dropwise until in excess.		
(iii) To the <b>third</b> part of the filtrate, add 2-3 drops barium nitrate solution followed by dilute nitric acid.		

<p>(iv) Use the <b>fourth</b> part of the filtrate to carryout a test of your choice to confirm the anion in <b>P</b>.</p> <p><b><u>TEST:</u></b></p>		
<p>(b) Wash the <b>residue</b> from part (a) with water thrice. Dissolve it in dilute hydrochloric acid. Add sodium hydroxide solution dropwise until in excess. Filter, keep both the <b>filtrate</b> and <b>residue</b>.</p>		
<p>(c) To the filtrate from part (b), add dilute hydrochloric acid until the solution is just acidic. Divide the resultant solution into <b>four</b> parts.</p>		
<p>(i) To the <b>first</b> part of acidic solution, add sodium hydroxide solution dropwise until in excess.</p>		
<p>(ii) To the <b>second</b> part of acidic solution, add 2-3 drops of potassium iodide solution.</p>		

(iii) To the <b>third</b> part of acidic solution, add ammonia solution dropwise until in excess		
(iv) Use the <b>fourth</b> part of acidic solution to carryout your own test to confirm the first cation in <b>P</b> . <b><u>TEST:</u></b>		
(d)Wash the <b>residue</b> from part (b) twice with sodium hydroxide solution. Dissolve it in dilute nitric acid. Divide the resultant solution into <b>four</b> parts.		
(i) To the <b>first</b> part of the solution, add sodium hydroxide solution dropwise until in excess.		
(ii) To the <b>second</b> part of the solution, add ammonia solution dropwise until in excess.		
(iii) To the <b>third</b> part of the solution, add one spatula endful of zinc dust and leave to stand.		

(iv) Use the <b>fourth</b> part of the solution to carryout a test of your choice to confirm the second cation in <b>P</b> .  <u><b>TEST:</b></u>		
---	--	--

(e i) Cations in **P**

are.....and.....

(ii) Anions in **P** are.....and.....

3. You are provided with an organic compound **Z**. You are required to carry out tests on **Z** to identify its nature. Identify any gas(es) evolved. Record your observations and deductions in **Table 4** below. (20 marks)

**Table 4**

TESTS	OBSERVATIONS	DEDUCTIONS
(a) Burn <b>Z</b> on a spatula end or porcelain dish.		
(b) Shake <b>two</b> spatula endfuls of <b>Z</b> with about 10cm <sup>3</sup> of water. Test the resultant solution with litmus paper. Divide the resultant solution into <b>six</b> parts.		



(i)	To the <b>first</b> part of the solution, add 2-3 drops of sodium carbonate solution.		
(ii)	To the <b>second</b> part of the solution, add 2-3 drops of Iron (III) chloride solution.		
(iii)	To the <b>third</b> part of the solution, add 2-3 drops of Brady's reagent.		
(iv)	To the <b>fourth</b> part of the solution, add 2-3 drops of acidified potassium manganate (VII) solution and heat.		
(v)	To the <b>fifth</b> part of the solution add Fehling's solution and heat.		
(vi)	To the <b>sixth</b> part of the solution, add 2-3 drops of Tollen's reagent and heat.		

(vii) To the <b>seventh</b> part of the solution, add 2-3 drops of Iodine solution followed by sodium hydroxide solution until the brown colour of Iodine is discharged. Warm the mixture and allow to cool.		
--	--	--

Comment on the nature of **Z**.

.....

.....

.....

.....

**END**

## **CONFIDENTIAL**

**FA1-** 0.1M HCl

**FA2-** 0.1M NaOH

**FA3-** 2M HCl

Solid **W-**  $\text{MgCO}_3$

**P-** Potassium chromate, Zinc carbonate and Copper (II) carbonate in ratio 1:3:2

**Z-** Glucose