

P525/3
Chemistry
Paper 3
July - August 2022
3 ¼ hours



UGANDA MUSLIM TEACHERS' ASSOCIATION

UMTA JOINT MOCK EXAMINATIONS 2022

Name

Centre/Index No. Signature

UGANDA ADVANCED CERTIFICATE OF EDUCATION

Chemistry

Paper 3

3 hours 15 Minutes

INSTRUCTIONS TO CANDIDATES:

- This paper consists of three compulsory questions.
- All questions must be answered in the spaces provided.
- Mathematical tables (3 – figure tables) and silent non-programmable scientific electronic calculators may be used.
- Candidates are advised to read through the paper and cross check with the apparatus and chemicals provided in the first fifteen minutes.

For Examiners' use only			
Q. 1	Q. 2	Q. 3	Total


- You are provided with the following;
FA1 which is 1M hydrochloric acid solution
FA2 which is approximately 1M sodium hydroxide solution
FA3 which is a 0.1M sulphuric acid
Solid T, which is impure tribasic acid H_3X

You are required to;

- Standardize solution FA2 using FA1
- Determine the percentage purity of solid T using FA2 and FA3.

Theory

Sodium hydroxide reacts with the acids according to the equations below:



Procedure A:

- Pipette $10cm^3$ of FA2 into a clean conical flask followed by 2 drops of phenolphthalein indicator and then titrate with FA1 from the burette until the end point. Repeat the titration until you obtain consistent readings. Enter your results in the table I below:

Table I

Volume of pipette used = 10.0 cm^3

Experiment	I	II	III
Final burette reading / cm^3	9.90	19.90	29.90
Initial burette reading / cm^3	0.00	10.00	20.00
Volume of FA 1 used / cm^3	9.90	9.90	9.90

(4 1/2 marks)

Titre values used to calculate average volume of FA1 (1/2 marks)

9.90 and 9.90

Therefore Average volume of FA1 = $\frac{9.90 + 9.90}{2} = 9.90$ cm^3 (2 1/2 marks)

Questions

Calculate the molar concentration of in FA2.

$1000cm^3$ of FA1 contains 1mole of HCl

$9.90cm^3$ of FA1 contains (1×9.90) moles of HCl = 9.9×10^{-3} moles

1mole of HCl reacted 1mole NaOH

9.9×10^{-3} mole of HCl reacted $(9.9 \times 10^{-3} \times 1)$ mole of NaOH = 9.9×10^{-3} mole

$10.0cm^3$ of FA2 contains 9.9×10^{-3} moles of NaOH

$1000cm^3$ of FA2 contains $(9.9 \times 10^{-3} \times 1000)$ moles of NaOH

= 0.99M

Procedure B

- b) Weigh accurately 3.4g of T and add about 50cm³ of water in a beaker. Stir to dissolve and transfer the contents of the beaker into a 250cm³ volumetric flask. Make up to the mark with distilled water. Label the resultant solution FA4.
- c) Pipette 10cm³ of FA4 into a conical flask. Add 10cm³ of FA2 using a measuring cylinder. Titrate the mixture with solution FA3 from the burette until the end point. Repeat the titration until you obtain consistent readings. Enter your results in table II below.

Table II

Mass of empty bottle + T	= 23.40 ✓		
Mass of empty bottle alone	= 20.0 ✓		
Mass of T alone	= 3.40 ✓		
Volume of pipette used =	10.0 ✓	cm ³	(2marks)
Experiment	I	II	III
Final burette reading /cm ³	39.10	49.00	44.00
Initial burette reading /cm ³	0.00	10.00	5.00
Volume of FA1 used /cm ³	39.10 ✓	39.00 ✓	39.00 ✓

Titre values used to calculate average volume of FA1 used are

39.00 and 39.00 cm³ (1/2 marks)

Therefore average volume of FA1 used = $\frac{39.00 + 39.00}{2}$ cm³

= 39.00 cm³ (2 1/2 marks)

Questions

d) Calculate the number of moles of;

i) sulphuric acid that reacted with the excess sodium hydroxide in FA4.

(1 1/2 marks)

1000cm³ of FA2 contains 0.1mole of H₂SO₄
 39.00cm³ of FA2 contains $\frac{0.1 \times 39.00}{1000}$ moles of H₂SO₄
 = 3.9×10^{-3} moles

(ii) Sodium hydroxide in 10cm³ of FA2 added to FA4 (1 1/2 marks)

$2\text{NaOH(aq)} + \text{H}_2\text{SO}_4\text{(aq)} \rightarrow \text{Na}_2\text{SO}_4\text{(aq)} + 2\text{H}_2\text{O(l)}$
 1mole of H₂SO₄ reacted with 2mole of NaOH
 3.9×10^{-3} moles of H₂SO₄ reacted with $2 \times 3.9 \times 10^{-3}$ moles of NaOH
 = 7.8 x 10⁻³ moles

ii) tribasic acid in 10cm^3 of FA4 that reacted with sodium hydroxide (1 mark)

$$\text{moles added} = (10 \times 0.99) \text{ moles} = 9.9 \times 10^{-3}$$

$$\text{moles of NaOH that reacted with T} = (9.9 \times 10^{-3} - 7.8 \times 10^{-3}) =$$

$$3 \text{ moles of NaOH reacted with 1 mole of T}$$

$$2.1 \times 10^{-3} \text{ moles of NaOH reacted with } (2.1 \times 10^{-3} \times \frac{1}{3}) \text{ moles} = 7.0 \times 10^{-4} \text{ moles}$$

e) Determine the percentage purity of T ($H=1$ $X=189$) (2½ marks)

$$100\text{cm}^3 \text{ of FA4 contains } 7.0 \times 10^{-4} \text{ moles of T}$$

$$250\text{cm}^3 \text{ of FA4 contains } (7.0 \times 10^{-4} \times 250) \text{ moles of T} = 0.0175 \text{ moles}$$

$$\text{Rfm of T} = (1 \times 3) + (189) = 192$$

$$1 \text{ mole of T weighs } 192\text{g}$$

$$0.0175 \text{ moles of T weigh } (192 \times 0.0175) \text{ g} = 3.36\text{g}$$

$$\% \text{age purity} = \frac{\text{mass of pure T}}{\text{mass of sample}} \times 100\%$$

$$= \frac{3.36}{3.4} \times 100\% = 98.8\%$$

2. You are provided with substance Y which contains three cations and one anion. You are required to carry out the following tests on Y to identify the cations and anion in it.

Identify any gases evolved. Record your observations and deductions in the table below.

Tests	Observations	Deductions
a). Heat a spatula end-ful of Y in a dry test tube until there is no further change.	Y is a mixture of white and brown crystals. - Colourless gas that turns moist paper blue and dense white fumes with HCl - Colourless gas that turns moist paper red and acidified $\text{K}_2\text{Cr}_2\text{O}_7$ from orange to green	Zn^{2+} , Al^{3+} , Fe^{3+} or Pb^{2+} NH_3 , NH_4^+ SO_2 , SO_4^{2-}
b) Shake two spatula end-ful of Y in a boiling tube with about 3cm^3 of water. Add dilute sodium hydroxide solution to the mixture drop wise until in excess. Warm and filter keep both the filtrate and the residue.	- Brown solution - Brown ppt insoluble in excess sodium hydroxide - Colourless gas that turns moist paper blue and dense white fumes with HCl - Brown residue - Colourless filtrate	Fe^{3+} Fe^{3+} NH_3 , NH_4^+ Fe^{3+} Zn^{2+} , Al^{3+} or Pb^{2+}

c) To the filtrate, add dilute nitric acid drop wise until the solution is just acidic. Divide the acidic solution into six parts. <i>Also dissolve around solution for soluble or dissolve.</i>	White precipitate soluble forming colorless solution	Zn^{2+}, Al^{3+} or Pb^{2+}
i) To the first part of the acidic solution, add dilute sodium hydroxide solution drop wise until in excess.	White precipitate soluble forming colorless solution	Zn^{2+}, Al^{3+} or Pb^{2+}
ii) To the second part of the acidic solution, add dilute ammonia solution drop-wise until in excess.	White precipitate insoluble	Al^{3+} or Pb^{2+}
iii) To the third part of the acidic solution, add 2-3 drops of potassium iodide solution.	No observable change	Pb^{2+} absent $\therefore Al^{3+}$ present <i>Allow Al^{3+} present.</i>
iv) To the fourth part of the acidic solution, add 2-3 drops of litmus solution followed by ammonia solution drop-wise until in excess. <i>Deny precipitate.</i>	Blue lake solution	Al^{3+} confirmed present
v) To the fifth part of the acidic solution, add 2-3 drops of lead (II) nitrate solution and heat.	White precipitate insoluble on heating	SO_4^{2-}
vi) Use the sixth part to carry out a test of your own choice to confirm the anion in Y Add $Ba(NO_3)_2$ followed by HNO_3	White precipitate insoluble in the acid.	SO_4^{2-} Confirmed
d) Wash the residue with water and dissolve it in dilute hydrochloric acid and divide the solution into three parts.	Brown solution	Fe^{3+}

All: $BaCl_2$ followed by HCl

with or without nitric acid

i) To the first part of the acidic solution, add dilute sodium hydroxide solution drop-wise until in excess.	Brown precipitate insoluble ✓	Fe^{3+} ✓
ii) To the second part of the acidic solution, add dilute ammonia solution drop wise until in excess.	Brown precipitate insoluble ✓	Fe^{3+} ✓
iii) To the third part of the acidic solution add 3-4 drops of potassium thiocyanate solution.	Blood-red coloration ✓	Fe^{3+} ✓ Confirmed

Allow ppt for precipitate.

The cations in Y are NH_4^+ ✓, Al^{3+} ✓ and Fe^{3+} ✓
The anion in Y is SO_4^{2-} ✓

3. You are provided with an organic compound Q. You are required to identify the nature of compound Q. Carry out the following tests on the compound and record your observations and deductions in the table below.

Tests	Observations	Deductions
2. • Burn a spatula end-ful of Q on a porcelain dish or at the end of a spatula.	Q burns with a blue yellow non-sooty flame ✓	Q is a saturated aliphatic organic compound of low carbon content ✓
3. • Shake 1cm^3 of Q with about 2cm^3 of water and test with litmus.	Miscible with water ✓ Has no effect on the litmus paper ✓	polar organic compound of low molecular mass ✓ Neutral organic compound ✓ Probably alcohol or aldehyde ✓
4. • To 0.5cm^3 of Q add 2-3 drops of sodium carbonate solution.	No effervescence or ✓ No bubbles ✓	Carboxylic acid absent ✓
2k. • To 0.5cm^3 of Q, add 2-3 drops of acidified potassium dichromate solution and heat.	The orange solution turns green ✓	Reducing agent ✓ Probably primary alcohol, secondary alcohol or aldehyde ✓

2	• To 0.5cm ³ of Q, add 2-3 drops of Brady's reagent.	No yellow precipitate or No observable change	Carbonyl compound absent.
2k	• To about 1cm ³ of Q, add acidified potassium dichromate solution and heat. Then add ethanol followed by 4-5 drops of concentrated sulphuric acid. Pour the mixture into a small beaker of cold water	Sweet & fruity smell. Ignore Sweet & fruity.	Ester formed Primary alcohol oxidised to Carboxylic acid
2	• To about 0.5cm ³ of Q, add about 4cm ³ of iodine solution followed by sodium hydroxide solution drop-wise until the brown color of iodine is just discharged. Warm the mixture and allow to stand.	Yellow precipitate	Primary al. Ethanol.
2	• To about 1cm ³ of Q, add about 5 drops of tollen reagent and heat gently	No Silver mirror ✓	Aldehyde absent

Comment on the nature of Q

1k Q is a saturated aliphatic primary alcohol with an methyl group on the same carbon as the hydroxyl group or Q is ethanol. ~~END~~***END***